


# 5 – 200 µl Filter tips

## 732612 | 732712



Dimensions in mm

<b>Volume</b>	5 – 200 µl
<b>Material</b>	PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	racked
<b>Quality</b>	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **



\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml