


0.1 – 20 µl Pipette tips

732122 | 732222



Dimensions in mm

Volume	0.1 – 20 µl
Material	PP, autoclavable at 121°C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml