


50 – 1000 µl Pipette tips

732112 | 732212



Dimensions in mm

| | |
|------------------|--|
| Volume | 50 – 1000 µl |
| Material | PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285 |
| Color | transparent |
| Packaging | racked |
| Quality | Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test ** |



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml