


5 – 200 µl Filter tips

732632 | 732732



Dimensions in mm

| | |
|------------------|--|
| Volume | 5 – 200 µl |
| Material | PP, Filter PE; autoclavable at 121°C (2 bar) acc. to DIN EN 285 |
| Color | transparent |
| Packaging | racked |
| Quality | Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP |



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml