



Product Data Sheet

LLG-Weighing boats 30 ml, PS white, diamond shape

FEATURES:	● Manufactured in high impact polystyrene ● Rounded corners provide for easy removal of powdered and granular samples. ● Recessed at the corners to facilitate controlled pour-off. ● Shelf life: 5 years from month of production on
STERILITY:	Non sterile
MATERIAL:	Polystyrene
THICKNESS:	0.18mm±10%/pc
WEIGHT:	1.1g±10%/pc
DIMENSIONS:	85x65mm
TREATMENT:	Antistatic

Lab Logistics Group GmbH Product
Management
Am Hambuch 1
D-53340 Meckenheim / Germany

Meckenheim, 2016-11-24

This document was sent online and is valid without signature

