




Product Data Sheet

LLG-Weighing boats, 100 ml



H (mm)	D (mm)	L (mm)	L1 (mm)
18.3	0.25	95.3	86.5

FEATURES: Manufactured in high impact Polystyrene and are ideal for handling both liquids and solids
 Rounded corners provide for easy removal of powdered and granular samples
 Recesses at the corners facilitate controlled pour-off

STERILITY: Non sterile

MATERIAL: Polystyrene

THICKNESS: 0.25mm±10%/pc

DIMENSIONS: 125mmx100mm

This document was sent online and is valid without signature