

TECHNICAL DATA SHEET
TC57333 SODIUM CHLORIDE pure
REV.07 - 16/07/2025

Product name	SODIUM CHLORIDE pure
Product code	TC57333
Formula – m.m. – C.A.S.	NaCl - 58,44 - 7647-14-5
TECHNICAL SPECIFICATIONS	
Appearance	White or colorless crystals
Sodium and Chloride Identification	Compliant
Assay	99,0 – 100,5 %
Alkalinity	≤ 0,50 ml HCl 0,01M
Bromides	≤ 100 ppm
Sulfates	≤ 200 ppm
Aluminum	≤ 0,2 ppm
Arsenic	≤ 1 ppm
Iron	≤ 2 ppm
Potassium	≤ 500 ppm
Heavy metals	≤ 3 ppm
Barium	≤ 5 ppm
Phosphates	≤ 25 ppm
Magnesium and alkali metals	≤ 100 ppm
Loss on drying	≤ 0,5 %
SUPPLIER SPECIFICATION DATA	
DOCUMENT ISSUED BY COMPUTER AND THEREFORE NOT SIGNED	