

TECHNICAL DATA SHEET	
TC86700 BUFFER solution pH 7,00 $\pm$ 0,02 at 25°C Colored YELLOW std	
REV. 01 DEL 14/06/2001	Control method MC384

Product Name	BUFFER solution pH 7,00 ± 0,02 at 25°C Colored YELLOW std	
Product Code	TC86700	
Composition	Solution of Potassium and Sodium Phosphate	
TECHNICAL SPECIFICATIONS		
Apparence	Clear yellow liquid	
pH (25° C)	6,98 – 7,02	
Correspondence Temperature (°C) / pH		
0	7,12	
5	7,07	
10	7,05	
15	7,03	
20	7,01	
25	7,00	
30	6,98	
35	6,97	
40	6,96	
45	6,96	
50	6,96	
The reported value is traceable through a chain of comparisons to the Standards Reference Material 186Ig-186Ilg of the National Institute of Standards and Technology (NIST - USA - www.nist.gov)		
DOCUMENT ISSUED BY COMPUTER AND THEREFORE NOT SIGNED		