

TECHNICAL DATA SHEET	
TC25700 IODINE 0.01N=N/100=0.005M	
REV. 4 - 20/08/2025	Control method MC136

Product name	IODINE 0.01N=N/100=0.005M
Product Code	TC25700
TECHNICAL SPECIFICATIONS	
Appearance	Clear red-brown liquid
Titration factor	0,995 - 1,005
<p>The value found is traceable through a chain of comparisons to the Standard Reference Material 136 of the National Institute of Standards and Technology (NIST – USA – www.nist.gov)</p>	
THIS DOCUMENT HAS BEEN MADE BY ELECTRONIC DATA PROCESSING AND THEREFORE NOT SIGNED	