

TECHNICAL DATA SHEET	
TC26000 IODINE 0,02N=N/50=0,01M	
REV. 04 – 08/10/2025	Control method MC138

Product name	IODINE 0,02N=N/50=0,01M
Product code	TC26000
TECHNICAL SPECIFICATIONS	
Appearance	Clear red-brown liquid
Titration factor	0,995 – 1,005
<p>The reported value is traceable through a chain of comparisons to the Reference Material product no. 42019 of the Federal Institute for Material Research and Testing (BAM - Germany - www.bam.de)</p>	
Expiry date	2 years from production in tightly closed bottle
THIS DOCUMENT HAS BEEN MADE BY ELECTRONIC DATA PROCESSING AND THEREFORE NOT SIGNED	