

TECHNICAL DATA SHEET	
TC34400 SODIUM HYDROXIDE 1N=N/1=1M	
REV. 06 – 16/10/2024	Control method MC204

Product name	SODIUM HYDROXIDE 1N=N/1=1M
Product code	TC34400
Composition	Reagent grade Sodium Hydroxide and reagenti water type III according to ASTM D1193
TECHNICAL SPECIFICATIONS	
Appearance	Clear colorless liquid
Titration factor	0,995 – 1,005
<p>The reported value is traceable through a chain of comparisons to the Reference Material product no. 93440 of the Federal Institute for Material Research and Testing (BAM - Germany - www.bam.de)</p>	
Expiry date	3 years from production in tightly closed bottle
THIS DOCUMENT HAS BEEN MADE BY ELECTRONIC DATA PROCESSING AND THEREFORE NOT SIGNED	