

<b>TECHNICAL DATA SHEET</b>	
TC48900 BROMOTHYMOL BLUE indicator 0.4% hydroalcoholic solution	
REV. 02 – 10/10/2025	Control method MC76

Product name	<b>BROMOTHYMOL BLUE indicator 0.4% hydroalcoholic solution</b>
Product code	<b>TC48900</b>
Composition	<b>Hydroalcoholic solution of Bromothymol Blue</b>
<b>TECHNICAL SPECIFICATIONS</b>	
<b>Appearance</b>	Clear dark green liquid
<b>Color change range (yellow-blue) pH</b>	5,8 – 7,6
<b>Sensitivity to use</b>	Conforme
<b>Suitability for use</b>	Acid/base indicator
<b>Expiry date</b>	5 years from production in tightly closed bottle
THIS DOCUMENT HAS BEEN MADE BY ELECTRONIC DATA PROCESSING AND THEREFORE NOT SIGNED	