

<b>TECHNICAL DATA SHEET</b>	
TC508I0 CARREZ I	
REV.3 - 29/04/2025	Control Method MC1236

Product name	<b>CARREZ I</b>	
Product code	<b>TC508I0</b>	
Composition	<b>Aqueous solution of Potassium Ferrocyanide</b>	
<b>TECHNICAL SPECIFICATIONS</b>		
<b>Appearance</b>	Clear yellowish liquid	
<b>Assay (Potassium Ferrocyanide)</b>	14,5 – 15,5 % w/v	
<b>Suitability for use:</b>	Oenology and food analysis in general	
<b>Expiry Date:</b>	5 years from production in tightly closed bottle	
DOCUMENT ISSUED BY COMPUTER AND THEREFORE NOT SIGNED		