

TECHNICAL DATA SHEET	
TC508II CARREZ II	
REV.2 - 29/04/2025	Control Method MC1236

Product name	CARREZ II	
Product code	TC508II	
Composition	Aqueous solution of Zinc Sulphate	
TECHNICAL SPECIFICATIONS		
Appearance	Clear, colorless liquid	
Assay (Zinc Sulphate Monohydrate)	29,5 – 30,5 % w/v	
Suitability for use:	Oenology and food analysis in general	
Expiry Date:	5 years from production in tightly closed bottle	
DOCUMENT ISSUED BY COMPUTER AND THEREFORE NOT SIGNED		