

TECHNICAL DATA SHEET	
TC53900 FEHLING B Potassium sodium tartrate alkaline	
REV. 04 – 10/10/2025	Control method MC404

Product name	FEHLING B Potassium sodium tartrate alkaline
Product code	TC53900
Composition	Aqueous solution of Potassium Sodium Tartrate and Sodium Hydroxide
TECHNICAL SPECIFICATIONS	
Appearance	Clear colorless liquid
Density (20°C)	1,26 – 1,28 g/ml
Alkalinity (NaOH)	119 – 121 g/l
Suitability for use	Sugar determination
Expiry date	3 years from production in tightly closed bottle
THIS DOCUMENT HAS BEEN MADE BY ELECTRONIC DATA PROCESSING AND THEREFORE NOT SIGNED	