

TECHNICAL DATA SHEET

TC54930 HYDROGEN PEROXIDE 30% very pure

REV.02 - 05/05/2025

Product name	HYDROGEN PEROXIDE 30% very pure
Product code	TC54930
Formula – m.m. – C.A.S.	H₂O₂ - 34,015 - 7722-84-1
TECHNICAL SPECIFICATIONS	
Aspect	Clear, colorless liquid
Assay	29,0 – 31,0 % p/p
Density	1,132 – 1,134 g/ml
Chlorides (Cl)	< 1 ppm
Phosphates (PO₄)	< 0,02%
Iron (Fe)	< 0,5 ppm
Lead (Pb)	< 0,1 ppm
Chrome (Cr)	< 0,5 ppm
Cobalt (Co)	< 0,1 ppm
Copper (Cu)	< 0,1 ppm
Zinc (Zn)	< 0,1 ppm
Reagent suitable for ASTM tests.	
SUPPLIER SPECIFICATION DATA	
DOCUMENT ISSUED BY COMPUTER AND THEREFORE NOT SIGNED	