

TECHNICAL DATA SHEET	
TC58711 INDICATOR for SULPHUR DIOXIDE in wines (OIV Method)	
REV. 02 – 10/10/2025	Control method MC2862

Product name	INDICATOR for SULPHUR DIOXIDE in wines - (OIV Method)	
Product code	TC58711	
Composition	Methylene Blue and Methyl Red in hydroalcoholic solution according to the OIV-MA-AS323-04A method	
TECHNICAL SPECIFICATIONS		
Appearance	Dark blue liquid	
Color change interval	5,0 – 5,5	
Suitability for use	Indicator for sulfur dioxide in wines	
Expiry date	4 years from production in tightly closed bottle	
THIS DOCUMENT HAS BEEN MADE BY ELECTRONIC DATA PROCESSING AND THEREFORE NOT SIGNED		