

Specification

**Albumin (BSA) Fraction V (pH 7.0) for Western blotting**

**A6588**

<b>Physical Description:</b>	Solid
<b>Product Code:</b>	A6588
<b>Product Name:</b>	Albumin (BSA) Fraction V (pH 7.0) for Western blotting
<b>Specifications:</b>	Proteases: not detectable Assay (Protein): min. 96 % pH (7 %; H <sub>2</sub> O; 20°C): 6.5 - 7.4 Heavy metals: max. 0.001 % NaCl: max. 2 % Immunoglobulins: not detectable Water (K.F.): max. 5 %
<b>WGK:</b>	1
<b>Storage:</b>	2 - 8°C
<b>Origin:</b>	from bovine serum
<b>M:</b>	approx. 68000 g/mol
<b>CAS:</b>	9048-46-8
<b>EINECS:</b>	232-936-2
<b>CS:</b>	35029070

**AppliChem GmbH**