

Grosseron being a part of you	FICHE TECHNIQUE		
	Date de création :	24/05/2023	
	Dernière Mise à jour :	14/10/2024	
	Version :	1	

Grosseron Collection	TUBE A CENTRIFUGER PP PINK ISAMPLE 50ML FD CONIQ STERILE 25/RACK BT500
Référence	GRR0501

DESCRIPTION

- Description 50ml Centrifuge Tubes, Bulk Pack, Sterilized
- Component Material/Resin PP/HDPE, free from heavy metal
- Colour Transparent ,Clear
- Dimension See drawing below
- Sterilization E-Beam SAL 10⁻⁶
- Quality Control Each lot is sampled and inspected in accordance with Standard Operating Procedures.
- Other Information For single use only


FEATURES

- Basic Features DNase/RNase Free, Non-pyrogenic
- Temperature Approx. -80 to +121 °C
- Autoclavability N/A
- Centrifugation, max, RCF 12000 x g
- Surface Treatment N/A
- Shelf Life 3 years after date of production

PACKAGING

- Package/Material PE
- Package Unit 25pcs/unit
- Units/Box 20 Unit/box
- Lot No. YYYY*** (year, serial number)

DRAWING



Ces spécifications ne peuvent être considérées que comme des informations générales et ne dispenseront pas l'utilisateur de son obligation de s'assurer que le produit est adapté pour l'utilisation prévue. En cas de doute, veuillez nous contacter.

FT GRR0501

Octobre 2024