



Référence

GRV0151

DESCRIPTION

- Description 15ml Centrifuge Tubes, Bulk Pack, Sterilized
- Component Material/Resin PP/HDPE, free from heavy metal
- Colour Transparent ,Clear
- Dimension See drawing below
- Sterilization E-Beam SAL 10⁻⁶
- Quality Control Each lot is sampled and inspected in accordance with Standard Operating Procedures.
- Other Information For single use only


FEATURES

- Basic Features DNase/RNase free, Non-pyrogenic
- Temperature Approx.-80 to+121°C
- Autoclavability N/A
- Centrifugation, max, RCF 12000 x g
- Surface Treatment N/A
- Shelf Life 3 years after date of production

PACKAGING

- Package/Material PE
- Package Unit 50pcs/unit
- Units/Box 10 Unit/box
- Lot No. YYYY*** (year, serial number)

DRAWING



Ces spécifications ne peuvent être considérées que comme des informations générales et ne dispenseront pas l'utilisateur de son obligation de s'assurer que le produit est adapté pour l'utilisation prévue. En cas de doute, veuillez nous contacter.