

	FICHE TECHNIQUE	
	Date de création :	27/04/2026
	Dernière Mise à jour :	27/04/2026
	Version :	1

	Pointe sans filtre universelle ISAMPLE 1000µL
Référence	GRV1000

DESCRIPTION	1000 ul Universal Fit Pipette Tips Bulk Bagged
- Component Material/Resin	PP(Polypropylene), 100% Virgin Medical Grade
- Colour	Blue
- Sterilization	N/A
- Quality Control	Each lot is sampled and inspected in accordance with Standard Operating Procedures.
- Other Information	For single use only

FEATURES	
- Basic Features	Non-pyrogenic Dnase/Rnase and human DNA Free,
- Temperature	Approx. -80°C to +121°C
- Autoclavability	Autoclavable at 121°C/15psi
- Shelf Life	5 years after date of production

PACKAGING	
- Package/Material	PE Bag, 1000/ Bag
- Package Unit	Bag/Case 1000 pcs/Bag , 5 Bags/5000 pcs/ Carton

Quality Control Criteria	
-Visual Observation:	No bending/Flash/Scratching/alien substance
-Dimensions	Compliance with drawings, within tolerance
-Packing	Quantity and Package Integrity
-Dnase/Rnase testing	Nuclease assay
-Compatibility	Eppendorf Pipettors, etc

Ces spécifications ne peuvent être considérées que comme des informations générales et ne dispenseront pas l'utilisateur de son obligation de s'assurer que le produit est adapté pour l'utilisation prévue. En cas de doute, veuillez nous contacter.

FT GRV1000

avril 26