



Acide trifluoroacétique

Synonyms:

TFA

CF₃COOH

Classification

Danger

Molecular weight

114.02 ONU

2699



CAS

76-05-1 Transport hazard class:

8

EEC-N:

200-929-3 Packing

I

Acide trifluoroacétique

> RS - Pour synthèse peptidique

Grade Spécifique

Identification (I.R.)	Conform	-
Density d _{20/4}	1.480 - 1.500	-
Colour	<= 10	APHA
Water (K.F.)	<= 200	ppm
Assay (acidimetric)	>= 99.7	% (w/w)
Chloride	<= 50	ppm



Sulfate	<= 50	ppm
Fluoride	<= 50	ppm
Identification (I.R.)	Conform	-
Density d20/4	1.480 - 1.500	-
Colour	<= 10	APHA
Water (K.F.)	<= 200	ppm
Assay (acidimetric)	>= 99.7	% (w/w)
Chloride	<= 50	ppm
Sulfate	<= 50	ppm
Fluoride	<= 50	ppm

ID	Size	Packaging
P0082103	100 ml	Flacon verre
P0082118	500 ml	Flacon verre
P0082147	2.5 l	Flacon verre sécurisé