



Ethanol absolu anhydre

Synonyms:

ETOH

C₂H₅OH**Classification****Danger**

Molecular weight

46 ONU**1170**

CAS

64-17-5 Transport hazard
class:**3**

EEC-N:

200-578-6 Packing**II****Ethanol absolu anhydre** > **RS - ERBA**dry - Solvants anhydres

Grade Specifique

Refractive index at 20°C	1.358 - 1.362	-
Water (K.F.)	<= 200	ppm
Colour	<= 10	APHA
Assay (alcoholic) at 20°C	>= 99.9	% (V/V)
Residue on evaporation	<= 10	ppm
Assay (GC)	>= 99.8	%
Acidity	<= 10	ppm



Aldehydes (as acetaldehyde)	<= 3	ppm
Esters (as CH3COOC2H5)	<= 25	ppm
Refractive index at 20°C	1.358 - 1.362	-
Water (K.F.)	<= 200	ppm
Colour	<= 10	APHA
Assay (alcoholic) at 20°C	>= 99.9	% (V/V)
Residue on evaporation	<= 10	ppm
Assay (GC)	>= 99.8	%
Acidity	<= 10	ppm
Aldehydes (as acetaldehyde)	<= 3	ppm
Esters (as CH3COOC2H5)	<= 25	ppm

ID	Size	Packaging
P013A10D03	100 ml	Flacon verre septum
P013A10D16	1 l	Flacon verre septum



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val de
Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System

