

Technical Specification

PRODUCT : Petroleum ether 40 - 60°C ERBAtec
CODE : P04599
CAS N° : 64742-49-0

TEST	U.M.	SPECIFICATION
Description	-	Coulourless liquid
Density at 20°C	-	0.64-0.67
Boiling point min.	°C	>40
Boiling point max.	°C	<70

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents w ould not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in