

50mL Centrifuge Tubes (Self-Standing)

Features

1. Made from virgin polypropylene to meet the general requirements of the laboratory.
2. Printed clear graduations with marking area
3. PP for the tube and HDPE for the cap
4. One-hand operation design
5. Sterilized by E-Beam
6. DNase, RNase and pyrogen free

Description

Temperature Range : Approx.-80 to+121 °C
 Shelf Life : 3 years after date of production
 Sterile : Yes

Product Range

Cat.No.	Volume (mL)	Packaging	Max.RCF (xg)	/Pack	/Case
602072	50	Bagged	5000	25	500

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.230, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA