

Product Name:	Cryo Boxe
Product Application:	For centralized placement of cryogenic vials.
Product Materials:	PC (cap) + ABS (bottom)
product Specifications:	Temperature range: stored at room temperature Shelf life: 5 Years after date of production (ensure package is in good)
Sterilization:	No
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: pass Packaging inspection: pass
Features	<ol style="list-style-type: none"> 1. The cap is made of polycarbonate plastic (PC), and the bottom is made of acrylonitrile butadiene styrene plastic (ABS). 2. Clear polycarbonate cover for easy identification. 3. Individually packed. 4. Alphanumeric grid for easy sample identification. 5. Each box includes drain holes and air vents to minimize condensation. 6. Withstand temperature range: Approx.-196°C to 121°C.


Product Range

Cat. No.	Product Description	/Pack	/Case
616224	For 0.5 ml vials	9*9	56
616325	For 1.2 ml vials	9*9	48
616021	For 2.0 ml vials	9*9	40
616421	For 3.0 ml vials	9*9	32
616121	For 4.0 ml vials	9*9	32
616521	For 5.0 ml vials	9*9	32
616204	For 0.5ml vials, Internal Thread Cryogenic Vial Only	10*10	56
616305	For 1.2ml vials, Internal Thread Cryogenic Vial Only	10*10	48
616001	For 2.0ml vials, Internal Thread Cryogenic Vial Only	10*10	40


Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Product Picture







Technical Drawing



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137




Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com


Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137


10*10 Cryo Boxes for 0.5 ml Vials
Internal Thread Cryogenic Vial Only



10*10 Cryo Boxes for 1.2 ml Vials
Internal Thread Cryogenic Vial Only



10*10 Cryo Boxes for 2.0 ml Vials
Internal Thread Cryogenic Vial Only



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)