

15mL Centrifuge Tubes

Features

1. Made from virgin polypropylene to meet the general requirements of the laboratory.
2. Printed clear graduations with marking area
3. PP for the tube and HDPE for the cap
4. One-hand operation design
5. Sterilized by E-Beam
6. DNase, RNase and pyrogen free

Description

Temperature Range : Approx. -80 to +121 °C
 Shelf Life : 3 years after date of production
 Sterile : Yes

Product Range

| Cat.No. | Volume (mL) | Packaging | Max.RCF (xg) | /Pack | /Case |
|---------|-------------|-----------|--------------|-------|-------|
| 601201 | 15 | PS Rack | 12000 | 50 | 500 |

Technical Drawing

