


250 mL & 500mL Centrifuge Tubes


The 250mL/500mL centrifuge tubes can process a large number of samples at one time and are suitable for separation and harvest of large cells, supernatants, bacteria, yeast and tissue samples, therefore they are widely used in the field of life sciences and clinical experiments and can meet the requirements of biological analysis.

Features

1. Unique structure design of the cap provides double security, further improve the seal ability.



2. Increased the thickness of conical bottom in order to improve the maximum centrifugal limit.



3. Polypropylene tube and high density polyethylene cap
4. Max. centrifugal limit 7000xg.
5. Temperature range: -80°C~120 °C
6. Clear scale
7. DNase, RNase and pyrogen free
8. Electron beam sterilization, SAL=10⁻⁶

Application


1. Separation of sediments after chemical reactions for harvest and application.
2. Harvest of biomacromolecules, inorganics and organics.
3. In biochemistry and other biological fields, usually used to centrifuge and collect a large volume of cells, biomacromolecules and supernatants.

NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.530, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA

Product Range

Cat. No.	Description	Unit
622001	250mL Centrifuge Tube, Bulk, Sterilized	6 /pk, 102 /cs
623001	500mL Centrifuge Tube, Bulk, Sterilized	6 /pk, 36 /cs
622002	250mL Centrifuge Tube, PS Rack, Sterilized	6 /pk, 24 /cs
623002	500mL Centrifuge Tube, PS Rack, Sterilized	4 /pk, 24 /cs

Technical Drawing



NEST China	Tel	: 0510-6878 8698 / 6878 8718	NEST USA	Tel	: +1-732-381-0268
	Email	: info@nest-wuxi.com		Email	: info@nestscientificusa.com
	Website	: www.cell-nest.com		Website	: www.nestscientificusa.com
	Add	: No.530, xida road, new district, Wuxi, Jiangsu, China		Add	: 1592 Hart St., Rahway, NJ07065, USA