

Lab Burner



Features

- ▶ Piezo ignition
- ▶ Max. temperature: 1300°C
- ▶ Adjustable flame
- ▶ Refillable butane tank
- ▶ International certification
CE certification
- ▶ One-Year Warranty

Applications

- ▶ Sterilizing
- ▶ Welding
- ▶ Heating
- ▶ Cooking

Ordering Information

- 177100-00**
Dragon 100, Torch Burner
- 177200-00**
Dragon 200, Lab Burner
- 177220-00**
Dragon 220, Lab Burner

	Dragon 100
	Dragon 200
	Dragon 220



Specification

Model	Dragon 100	Dragon 200	Dragon 220
Max. temperature	1300°C	1300°C	1300°C
Fuel tank capacity	40 g	35 g	35 g
Working time	60 min	60 min	60 min
Safe switch	Yes	Yes	Yes
Dimension (LxWxH)	7.9 x 13.7 x 16.1 cm	7.9 x 11.5 x 9.2 cm	7.8 x 11.5 x 9.4 cm
Net weight	380 g	407 g	425 g